








Winston-Lutz Analysis


Winston-Lutz Analysis

=====


Number of images: 12
Maximum 2D CAX->BB distance: 0.67mm
Median 2D CAX->BB distance: 0.31mm
Shift to iso: facing gantry, move BB: RIGHT 0.10mm; IN 0.25mm; UP 0.04mm
Gantry 3D isocenter diameter: 0.61mm (4/12 images considered)
Maximum Gantry RMS deviation (mm): 0.17mm
Maximum EPID RMS deviation (mm): 0.68mm
Gantry+Collimator 3D isocenter diameter: 0.69mm (8/12 images considered)
Collimator 2D isocenter diameter: 0.28mm (5/12 images considered)
Maximum Collimator RMS deviation (mm): 0.52
Couch 2D isocenter diameter: 0.62mm (5/12 images considered)
Maximum Couch RMS deviation (mm): 0.67



Gantry images



Collimator images



Couch images

