








Winston-Lutz Analysis





Winston-Lutz Analysis

Number of images: 12
 Maximum 2D CAX->BB distance: 0.67mm
 Median 2D CAX->BB distance: 0.31mm
 Shift to iso: facing gantry, move BB: RIGHT 0.10mm; IN 0.25mm; UP 0.04mm
 Gantry 3D isocenter diameter: 0.61mm (4/12 images considered)
 Maximum Gantry RMS deviation (mm): 0.17mm
 Maximum EPID RMS deviation (mm): 0.68mm
 Gantry+Collimator 3D isocenter diameter: 0.69mm (8/12 images considered)
 Collimator 2D isocenter diameter: 0.28mm (5/12 images considered)
 Maximum Collimator RMS deviation (mm): 0.52
 Couch 2D isocenter diameter: 0.62mm (5/12 images considered)
 Maximum Couch RMS deviation (mm): 0.67




Winston-Lutz Analysis

Gantry images



Winston-Lutz Analysis

Collimator images



Winston-Lutz Analysis

Couch images

